

## Anti-Interferon Alpha 1 antibody [2-48] (STJ16103684)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA/IHC-F/WB     |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                |
|---------------------|------------------------------------------------|
| Clonality           | Monoclonal                                     |
| Concentration       | 100 Mu g/mL                                    |
| Conjugation         | Unconjugated                                   |
| Formulation         | PBS with 0.02% sodium azide                    |
| Isotype             | IgG1                                           |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|