

Anti-IgG3 antibody [LO-MG3-13] (STJ16103688)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/ELISA
Host / Source	Rat
Reactivity	Mouse/Swine/Horse/Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	PBS
Isotype	IgMk
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---