

Anti-Double stranded DNA antibody [AE-2] (STJ16103689)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC-F/IHC-P/IF
Host / Source	Mouse
Reactivity	Human/Mouse/Rat/Rabbit/All Eukaryotes

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Dilution Range	FC 1-2ug million cells IHC-F IHC-P IF all 1-2 ug ml
Formulation	PBS with 0.02% sodium azide
Isotype	IgG3
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---