

Anti-MUC5AC antibody [45M1] (STJ16103701)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IHC-F/IHC-P/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat/Rabbit/Monkey/Chicken/Cow/Cat/Pig/Hedgehog

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	4586
Gene Symbol	MUC5AC
UniProt ID	MUC5A_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---