

Anti-RBP antibody [G4E4] (STJ16103717)

GENERAL INFORMATION

| | |
|---------------|------------------------------------------------------|
| Product Type | Primary antibodies |
| Applications | ELISA/IHC-F/IHC-P/WB |
| Host / Source | Mouse |
| Reactivity | Chimpanzee/Goat/Human/Monkey/Mouse/Rabbit/Rat/Rhesus |

PRODUCT PROPERTIES

| | |
|---------------------|------------------------------------------------|
| Clonality | Monoclonal |
| Concentration | 100 Mu g/mL |
| Conjugation | Unconjugated |
| Formulation | PBS with 0.02% sodium azide |
| Isotype | IgG1 |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

ADDITIONAL INFORMATION

| | |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|