

## Anti-HLA-B27 antibody [EP-4] (STJ16103721)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/FC/IF
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide
Isotype	IgMk
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---