

Anti-Fibronectin antibody [616] (STJ16103740)

GENERAL INFORMATION

| | |
|---------------|----------------------|
| Product Type | Primary antibodies |
| Applications | FC/IHC-F/IHC-P/IF/WB |
| Host / Source | Mouse |
| Reactivity | Human/Dog/Ferret/Pig |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 100 Mu g/mL |
| Conjugation | Unconjugated |
| Formulation | PBS with 0.02% sodium azide |
| Isotype | IgG2ak |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 2335 |
| Gene Symbol | FN1 |
| UniProt ID | FINC_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|