

## Anti-Myeloid Specific Antigen antibody [BM-2] (STJ16103742)

### GENERAL INFORMATION

|               |                      |
|---------------|----------------------|
| Product Type  | Primary antibodies   |
| Applications  | FC/IHC-F/IHC-P       |
| Host / Source | Mouse                |
| Reactivity    | Human/Macaque Monkey |

### PRODUCT PROPERTIES

|                     |                                                |
|---------------------|------------------------------------------------|
| Clonality           | Monoclonal                                     |
| Concentration       | 100 Mu g/mL                                    |
| Conjugation         | Unconjugated                                   |
| Formulation         | PBS with 0.02% sodium azide                    |
| Isotype             | IgG1k                                          |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|