

Anti-Laminin γ antibody [A5] (STJ16103749)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC-F/IHC-P/IF
Host / Source	Rat
Reactivity	Human/Monkey/Mouse/Pig/Rat/Rhesus

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 μ g/mL
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide
Isotype	IgG2ak
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---