

## Anti-CD123 antibody (26-171) [ZR425] (STJ180630)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-P              |
| Host / Source | Rabbit             |
| Reactivity    | Human              |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                                                                                                                    |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                                                                                                                                         |
| Conjugation         | Unconjugated                                                                                                                                                                                                                                       |
| Purification        | Affinity purified                                                                                                                                                                                                                                  |
| Dilution Range      | 1:100-200                                                                                                                                                                                                                                          |
| Formulation         | Tris-HCl buffer containing stabilizing protein (BSA) and <0.1% ProClin                                                                                                                                                                             |
| Isotype             | IgG                                                                                                                                                                                                                                                |
| Storage Instruction | Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C. |

### TARGET INFORMATION

|                  |                                        |
|------------------|----------------------------------------|
| Gene ID          | <a href="#">3563</a>                   |
| Gene Symbol      | <a href="#">IL3RA</a>                  |
| UniProt ID       | <a href="#">IL3RA_HUMAN</a>            |
| Immunogen Region | 26-171                                 |
| Specificity      | Positive Control: Lymph node or thymus |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti-CD123 antibody IHC Human tonsil  
 STJ180630