

Anti-STMN1 antibody (1-100 aa) [ABT347] (STJ190092)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/WB
Host / Source	Mouse
Reactivity	Human/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution Range	IHC-P 1:200-400 WB 500-1000
Formulation	Liquid in PBS pH7.2, 0.03% Proclin 300, with stabilizing protein.
Isotype	IgG2ak
Molecular Weight	Calculated: 17kD Observed: 17kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	3925
Gene Symbol	STMN1
UniProt ID	STMN1_HUMAN
Immunogen Region	1-100 aa
Specificity	The antibody can specifically recognize human Stathmin protein. In western blotting of Hela cell lysate, the antibody can label a 17 kDa band corresponding to Stathmin.. The antibody was also Predict

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Human tonsil tissue was stained with Anti-Stathmin (ABT347) Antibody



Human tonsil tissue was stained with Anti-Stathmin (ABT347) Antibody



Various whole cell lysates were separated by 10% SDS-PAGE, and the membrane was blotted with anti-Stathmin (ABT347) antibody. The HRP-conjugated Goat anti-mouse IgG (H + L) antibody was used to detect the antibody. Lane 1: Neuro-2a Lane 2: Hela Lane 3: Jurkat Lane 4: NIH-3T3 Predicted band size: 17kDa Observed band size: 17kDa