

## Anti-SLC6A18 antibody (STJ191430)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/ELISA           |
| Host / Source | Rabbit             |
| Reactivity    | Human/Rat/Mouse    |

### PRODUCT PROPERTIES

|                     |                                                                                                                       |
|---------------------|-----------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                            |
| Concentration       | 1 mg/mL                                                                                                               |
| Conjugation         | Unconjugated                                                                                                          |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Dilution Range      | WB 1:500-2000<br>ELISA 1:5000-20000                                                                                   |
| Formulation         | Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.                                                         |
| Isotype             | IgG                                                                                                                   |
| Molecular Weight    | Observed: 69kD                                                                                                        |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                        |

### TARGET INFORMATION

|             |                                                                 |
|-------------|-----------------------------------------------------------------|
| Gene ID     | <a href="#">348932</a>                                          |
| Gene Symbol | <a href="#">SLC6A18</a>                                         |
| UniProt ID  | <a href="#">S6A18_HUMAN</a>                                     |
| Specificity | S6A18 Polyclonal Antibody detects endogenous levels of protein. |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|