

## Anti-NR1H4 antibody (STJ192265)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/ELISA           |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES


|                     |                                                                                                                       |
|---------------------|-----------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                            |
| Concentration       | 1 mg/mL                                                                                                               |
| Conjugation         | Unconjugated                                                                                                          |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Dilution Range      | WB 1:500-2000<br>ELISA 1:5000-20000                                                                                   |
| Formulation         | Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.                                                         |
| Isotype             | IgG                                                                                                                   |
| Molecular Weight    | Observed: 53kD                                                                                                        |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                        |

### TARGET INFORMATION

|             |                                                                 |
|-------------|-----------------------------------------------------------------|
| Gene ID     | <a href="#">9971</a>                                            |
| Gene Symbol | <a href="#">NR1H4</a>                                           |
| UniProt ID  | <a href="#">NR1H4_HUMAN</a>                                     |
| Specificity | NR1H4 Polyclonal Antibody detects endogenous levels of protein. |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Huang, Jianbo, et al. "Effects of mulberry leaf on experimental hyperlipidemia rats induced by high-fat diet. " Experimental and therapeutic medicine 16. 2 (2018) : 547-556.