

Anti-HTR3B antibody (140-220 aa) (STJ192609)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Rat/Mouse
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	1 mg/mL
	Conjugation	Unconjugated
	Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
	Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
	Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
	Isotype	IgG
	Molecular Weight	Observed: 48kD
	Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	9177
	Gene Symbol	HTR3B
	UniProt ID	5HT3B_HUMAN
	Immunogen Region	140-220 aa
	Specificity	5HT3B Polyclonal Antibody detects endogenous levels of protein.
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.