

## Anti-LRRC3B antibody (115-165 aa) (STJ194951)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC             |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse        |

### PRODUCT PROPERTIES


|                     |                                                                                                                       |
|---------------------|-----------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                            |
| Concentration       | 1 mg/mL                                                                                                               |
| Conjugation         | Unconjugated                                                                                                          |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Dilution Range      | WB 1:500-2000<br>IHC-P 1:50-300                                                                                       |
| Formulation         | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.                                               |
| Isotype             | IgG                                                                                                                   |
| Molecular Weight    | Calculated: 28kD                                                                                                      |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                        |

### TARGET INFORMATION


|                  |                                                                 |
|------------------|-----------------------------------------------------------------|
| Gene ID          | <a href="#">116135</a>                                          |
| Gene Symbol      | <a href="#">LRRC3B</a>                                          |
| UniProt ID       | <a href="#">LRC3B_HUMAN</a>                                     |
| Immunogen Region | 115-165 aa                                                      |
| Specificity      | This antibody detects endogenous levels of LRC3B at Human/Mouse |

### ADDITIONAL INFORMATION

**Note** STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of lysates from THP-1 cells, primary antibody was diluted at 1:1000, 4°C overnight STJ194951



Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200 (4°C overnight). STJ194951