

Anti-SORBS2 antibody (794-844 aa) (STJ195823)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Calculated: 121kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	8470
Gene Symbol	SORBS2
UniProt ID	SRBS2_HUMAN
Immunogen Region	794-844 aa
Specificity	This antibody detects endogenous levels of SRBS2 at Human/Mouse/Rat

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



miR-21-3p in extracellular vesicles from vascular fibroblasts of spontaneously hypertensive rats promotes proliferation STJ195823



Western blot analysis of lysates from HUVECs, primary antibody was diluted at 1:1000, 4° overnight STJ195823