

Anti-PPARG antibody [Mix] (STJ196657)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution Range	IHC-P 1:50-200 WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Molecular Weight	Observed: 53-57kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	5468
Gene Symbol	PPARG
UniProt ID	PPARG_HUMAN
Specificity	This antibody detects endogenous levels of PPAR Gamma at Human, Mouse, Rat

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunohistochemical analysis of paraffin-embedded Human Placenta Tissue using PPAR γ mouse Monoclonal antibody diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Placenta Tissue using PPAR γ mouse Monoclonal antibody diluted at 1:200.



Western blot analysis of 1) HeLa Cell, 2) 293T Cell lysate using PPAR γ mouse Monoclonal antibody diluted at 1:2,000.