

Anti-PTH antibody (50-115 aa) [PT0210] (STJ196888)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/IF
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution Range	IHC-P 1:100-500 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG1k
Molecular Weight	Calculated: 13kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	5741
Gene Symbol	PTH
UniProt ID	PTHY_HUMAN
Immunogen Region	50-115 aa
Specificity	This antibody detects endogenous levels of human Parathyroid Hormone (PTH). Heat-induced epitope retrieval (HIER) TRIS-EDTA of pH8.0 was highly recommended as antigen repair method in paraffin section

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Human parathyroidoma tissue was stained with Anti-Parathyroid Hormone (ABT215) Antibody



Human parathyroidoma tissue was stained with Anti-Parathyroid Hormone (ABT215) Antibody



Human parathyroidoma tissue was stained with Anti-Parathyroid Hormone (ABT215) Antibody



Human parathyroidoma tissue was stained with Anti-Parathyroid Hormone (ABT215) Antibody