

Anti-CD14 antibody (1-100 aa) [ABT021] (STJ197096)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/WB/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution Range	IHC-P 1:100-500 WB 1:200-1000 IF 1:100-500
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG2bk
Molecular Weight	Calculated: 40kD Observed: 40kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	929
Gene Symbol	CD14
UniProt ID	CD14_HUMAN
Immunogen Region	1-100 aa
Specificity	The antibody can specifically recognize human CD14 protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Human adrenal gland tissue was stained with anti-CD14 (ABT021) antibody.



Human kidney tissue was stained with anti-CD14 (ABT021) antibody.



Human tonsil tissue was stained with anti-CD14 (ABT021) antibody.