

Anti-GYPA antibody (1-100 aa) [ABT179] (STJ197158)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution Range	IHC-P 1:100-500 IF 1:100-500
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG1k
Molecular Weight	Calculated: 16kD Observed: 37kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	2993
Gene Symbol	GYPA
UniProt ID	GLPA_HUMAN
Immunogen Region	1-100 aa
Specificity	The antibody can specifically recognize human Glycophorin A protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Human acute myeloid leukemia tissue was stained with Anti-Glycophorin A (ABT179) Antibody



Human kidney tissue was stained with Anti-Glycophorin A (ABT179) Antibody



Human placenta tissue was stained with Anti-Glycophorin A (ABT179) Antibody



Human spleen tissue was stained with Anti-Glycophorin A (ABT179) Antibody