

## Anti-ACP5 antibody (250-325 aa) [ABT522] (STJ197195)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/IF
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution Range	IHC-P 1:100-500 IF 1:100-500
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG1k
Molecular Weight	Calculated: 37kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION


Gene ID	<a href="#">54</a>
Gene Symbol	<a href="#">ACP5</a>
UniProt ID	<a href="#">PPA5_HUMAN</a>
Immunogen Region	250-325 aa
Specificity	The antibody can specifically recognize human TRAP protein.. The antibody was also Predict react with Mouse;Rat

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Human giant cell tumor tissue was stained with Anti-TRAP (ABT522) Antibody



Human lung tissue was stained with Anti-TRAP (ABT522) Antibody