

Anti-TYR antibody (250-350 aa) [ABT242] (STJ197248)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution Range	IHC-P 1:200-400
Formulation	Liquid in PBS pH7.2, 0.03% Proclin 300, with stabilizing protein.
Isotype	IgG1k
Molecular Weight	Calculated: 60kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	7299
Gene Symbol	TYR
UniProt ID	TYRO_HUMAN
Immunogen Region	250-350 aa
Specificity	The antibody can specifically recognize human Tyrosinase protein.. The antibody was also Predict react with Mouse;Rat

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Human malignant melanoma tissue was stained with Anti-Tyrosinase (ABT242) Antibody



Human malignant melanoma tissue was stained with Anti-Tyrosinase (ABT242) Antibody



Human malignant melanoma tissue was stained with Anti-Tyrosinase (ABT242) Antibody



Human skin tissue was stained with Anti-Tyrosinase (ABT242) Antibody