

Anti-KCNIP3 antibody (169-219 aa) (STJ197369)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/IF
Host / Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	IHC 1:100-200 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 29kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	30818
Gene Symbol	KCNIP3
UniProt ID	CSEN_HUMAN
Immunogen Region	169-219 aa
Specificity	KChIP3 protein (A250) detects endogenous levels of KChIP3

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunohistochemical analysis of paraffin-embedded Rat BrainTissue using KChIP3 Rabbit pAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded mouse BrainTissue using KChIP3 Rabbit pAb diluted at 1:200.