

Anti-AKT3 antibody (STJ197445)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS pH7.4, 0.5% BSA, 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG
Molecular Weight	Observed: 56kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10000
Gene Symbol	AKT3
UniProt ID	AKT3_HUMAN
Specificity	The antibody detects endogenous Akt3 protein.

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



The picture was kindly provided by our customer. primary antibody was diluted at 1:5000. Loading control antibody was diluted at 1:20000

Western blot analysis of 1) HeLa, 2) 3T3, 3) Mouse Brain, 4) Rat Brain tissue, diluted at 1:2000. Secondary antibody (STJ99666) was diluted at 1:20000