

Anti-RAET1G antibody (192-242 aa) (STJ197610)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB
	Host / Source	Rabbit
	Reactivity	Human/Rat/Mouse
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	1 mg/mL
	Conjugation	Unconjugated
	Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
	Dilution Range	WB 1:1000-2000
	Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
	Isotype	IgG
	Molecular Weight	Calculated: 37kD
	Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	353091
	Gene Symbol	RAET1G
	UniProt ID	ULBP5_HUMAN
	Immunogen Region	192-242 aa
	Specificity	This antibody detects endogenous levels of human RET1G
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.