

Anti-Phospho-HDAC8-S39 antibody (STJ22168)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/ICC
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF/ICC:1:100-1:200
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 42kDa Observed Mw:
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	55869
Gene Symbol	HDAC8
UniProt ID	HDAC8_HUMAN
Specificity	A phospho specific peptide corresponding to residues surrounding S39 of human HDAC8

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunohistochemistry of paraffin-embedded human lung carcinoma using Phospho-HDAC8-S39 antibody (STJ22168).



Immunofluorescence analysis of methanol-fixed HeLa cells using Phospho-HDAC8-S39 antibody (STJ22168).