

Anti-Phospho-IGF1R/INSR-Y1150/Y1135 antibody (STJ22192)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA: Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 154kDa/155kDa/156kDa Observed Mw: 80kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3643 3480
Gene Symbol	INSR IGF1R
UniProt ID	INSR_HUMAN IGF1R_HUMAN
Immunogen Sequence	TDYYR
Specificity	A synthetic phosphorylated peptide around Y1150 of human IGF1R/INSR (NP_001073285.1).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of lysates from 293 and HepG2 cells, using Phospho-INSR-Y1150/IGF1R-Y1135 Rabbit pAb (STJ22192) at 1:1000 dilution. 293 cells were treated by EGF (25 Mu g/mL) for 30 minutes after serum-starvation overnight. HepG2 cells were treated by EGF (100ng/mL) for 30 minutes after serum-starvation overnight. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% BSA.