

Anti-ABCF2 antibody (1-250) (STJ22464)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IF/CC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 71kDa Observed Mw: 72kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	10061
Gene Symbol	ABCF2
UniProt ID	ABCF2_HUMAN
Immunogen Region	1-250
Immunogen Sequence	MPSDLAKKKAAKKKEAAKAR QRPRKGHEENGDDVVTEPQVA EKNEANGRETTEVDLLTKEL EDFEMKKAARAVTGVLASH PNSTDVHIIINLSLTFHGQEL LSDTKLELNSGRRYGLIGLN GIGKSMLLSAIGKREVPPIE HIDIYHLTREMPSPDKTPLH CVMEVDTERAMLEKEAERLA HEDAECEKLMELYERLEELD ADKAEMRASRLHGLGFTPA MQRKKLKDFSGGWRMRVAL
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human ABCF2 (NP_005683.2).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunofluorescence analysis of U2OS cells using ABCF2 antibody (STJ22464) at dilution of 1:100. Blue: DAPI for nuclear staining.



Western blot analysis of extracts of various cell lines, using ABCF2 antibody (STJ22464) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.