

Anti-AKR7A2 antibody (100-359) (STJ22570)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 40kDa Observed Mw: 35kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8574
Gene Symbol	AKR7A2
UniProt ID	ARK72_HUMAN
Immunogen Region	100-359
Immunogen Sequence	VKIATKANPWDGKSLKPDSV RSQLETSLKRLQCPQVDFLY LHAPDHGTPVEETLHACQRL HQEGKFVELGLSNYASWEVA EICTLCKSNGWILPTVYQGM YNATTRQVETELFPCLRHFGL RLFYAYNPLAGLLTGKYKY EDKDGKQPVGRFFGNSWAET YRNRFWKEHHFEIALVEKA LQAAYGASAPSVTSAALRWM YHHSQLQGAHGDAVILGMSS LEQLEQNLAATEEGPLEPA
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 100-359 of human AKR7A2 (NP_003680.2).

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Western blot analysis of extracts of various cell lines, using AKR7A2 antibody (STJ22570) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit.