

Anti-ATG13 antibody (281-517) (STJ22706)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 IF/CC:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 57kDa Observed Mw: 72kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	9776
Gene Symbol	ATG13
UniProt ID	ATG13_HUMAN
Immunogen Region	281-517
Immunogen Sequence	ADLAYPVVFAAGLNATHPHQ LMVPGKEGGVPLAPNQPVHG TQADQERLATCTPSDRTHCA ATPSSSEDTETYSNSSEGRA SPHDVLETIFVRKVGAFVVK PINQVTLTSLDIPFAMFAPK NLELEDTPMVNPPDSPETE SPLQGLSHSDGSSGGSSGNT HDDFVIMDFKPAFSKDDILP MDLGTFYREFQNPQLSSLS
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 281-480 of human ATG13 (NP_001136145.1).

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of various lysates, using ATG13 antibody (STJ22706) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% non-fat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 25s.



Immunofluorescence analysis of HeLa cells using ATG13 rabbit polyclonal antibody (STJ22706) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Immunofluorescence analysis of NIH/3T3 cells using ATG13 rabbit polyclonal antibody (STJ22706) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.