

Anti-DTNBP1 antibody (1-351) (STJ23439)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 39kDa Observed Mw: 48kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84062
Gene Symbol	DTNBP1
UniProt ID	DTBP1_HUMAN
Immunogen Region	1-351
Immunogen Sequence	MLETLRERLLSVQQDFTSGL KTLSDKSREAKVKSKPRTVP FLPKYSAGLELLSRYEDTWA ALHRRAKDCASAGELVDSEV VMLSAHWKKTSLVLEQEQ LQQLPALIADLESMTANLTH LEASFEEVENLLHLEDLCG QCELERCKHMQSQLENYKK NKRKELETFFKAELDAEHAQK VLEMEHTQMKLKERQKFFE EAFQQDMEQYLSTGYLQIAE RREPIGSMSSMEVNVDMLE
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-351 of human DTNBP1 (NP_115498.2).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of extracts of various cell lines, using DTNBP1 antibody (STJ23439) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST.