

Anti-EIF2S1 antibody (1-260) (STJ23509)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:5000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 36kDa Observed Mw: 36kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1965
Gene Symbol	EIF2S1
UniProt ID	IF2A_HUMAN
Immunogen Region	1-260
Immunogen Sequence	MPGLSCRFYQHKFVEDEVV MVNVRSIAEMGAYVSLLEYN NIEGMILLSELSRRRRIRSIN KLIRIGRNECVVIRVDKEK GYIDLSKRRVSPPEAIKCED KFTKSKTVYSILRHVAEVL YTKDEQLLESLFQRTAWVFDD KYKRPGYGAYDAFKHAVSDP SILDSLNLNEDEREVLINNI NRRLTPQAVKIRADIEVACY GYEGIDAVKEALRAGLNCST ENMPIKINLIAPPRYVMTT
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human eIF2 Alpha (NP_004085.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of various lysates, using eIF2 Alpha antibody (STJ23509) at 1:400 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 15s.