

Anti-ENO3 antibody (1-190) (STJ23540)

GENERAL INFORMATION

	Product Type	Primary antibodies
	Applications	WB/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:4000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 47kDa Observed Mw: 47kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2027
Gene Symbol	ENO3
UniProt ID	ENOB_HUMAN
Immunogen Region	1-190
Immunogen Sequence	MAMQKIFAREILDSRGNPTV EVDLHTAKGRFRAAVPSGAS TGIYEALELRDGDKGRYLGK GVLKAVENINNTLGPALLQK KLSVVDQEKVDKFMIELDGT ENKSKFGANAILGVSLAVCK AGAAEKGVPLYRHIADLAGN PDLILPVPAPFNVINGGSHAG NKLAMQEFMILPVGASSFKE AMRIGAEVYH
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human ENO3 (NP_001967.3).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of extracts of various cell lines, using ENO3 antibody (STJ23540) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.