

Anti-ENPP2 antibody (564-863) (STJ23543)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 99kDa Observed Mw: 110-115kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5168
Gene Symbol	ENPP2
UniProt ID	ENPP2_HUMAN
Immunogen Region	564-863
Immunogen Sequence	DLGCTCDDKVEPKNKLDELN KRLHTKGSTEERHLLYGRPA VLYRTRYDILYHTDFESGYS EIFLMPLWTSYTVSKQAEVS SVPDHLTSCVRPDVVRVSPSF SQNCLAYKNDKQMSYGFLFP PYLSSSSPEAKYDAFLVTNMV PMYPAFKRVWNYFQRLVKK YASERNGVNVISGPIFDYDY DGLHDTEDKIKQYVEGSSIP VPTHYYSIITSCLDFTQPAD KCDGPLSVSSFILPHRPDN
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 564-863 of human ENPP2 (NP_001035181.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of extracts of various cell lines, using ENPP2 antibody (STJ23543) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.