

Anti-FAAH antibody (450-579) (STJ23602)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 63kDa Observed Mw: 60kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2166
Gene Symbol	FAAH
UniProt ID	FAAH1_HUMAN
Immunogen Region	450-579
Immunogen Sequence	EIEVYRKTIVIAQWRALDLDV VLTPLAPALDLNAPGRATG AVSYTMLYNCLDFPAGVVPV TTVTAEDAEQMEHYRGYFGD IWDKMLQKGMKKSGLPVAV QCVALPWQEELCLRFMREVE RLMTPEKQSS
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 450-579 of human FAAH11 (NP_001432.2).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of extracts of various cell lines, using FAAH11 antibody (STJ23602) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST.