

Anti-GNAS antibody (1-100) (STJ23816)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 28kDa/44kDa/45kDa/77kDa/109kDa/111kDa Observed Mw: 43kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2778 2778 2778 GNAS GNAS GNAS1_HUMAN GNAS3_HUMAN GNAS2_HUMAN
Immunogen Region	1-100
Immunogen Sequence	MGCLGNSKTEDQRNEEKAQR EANKKIEKQLQKDKQVYRAT HRLLLLGAGESGKSTIVKQM RILHVNGFN GEGGEEDPQAA RSNSDGEKATKVQDIKNNLK
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GNAS (P63092).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of extracts of various cell lines, using GNAS antibody (STJ23816) at 1:400 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.