

Anti-MAD1L1 antibody (360-700) (STJ24440)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IF/CC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 83kDa Observed Mw: 95kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	8379
Gene Symbol	MAD1L1
UniProt ID	MD1L1_HUMAN
Immunogen Region	360-700
Immunogen Sequence	RGLEKARQQLQEELRQVSGQ LLEERKKRETHEALARRLOK RVLLLTKERDGMRAILGSYD SELTPAEYSPQLTRRMREAE DMVQKVHSHSAEMEQLSQA LEELGGQKQRADMLEMELKM LKSQSSSAEQSFLFSREEAD TLRKVEELEGERSRLEEEK RMLEAQLERRALQGDYDQSR TKVLHMSLNPTSVARQRLRE DHSQLQAECERLRGLLRAME RGGTVPADLEAAAAASLPSS
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 360-700 of human MAD1/MAD1L1 (NP_001013859.1).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunofluorescence analysis of U2OS cells using MAD1/MAD1L1 antibody (STJ24440).



Western blot analysis of extracts of various cell lines, using MAD1/MAD1L1 antibody (STJ24440) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit.