

Anti-NFATC1 antibody (20-300) (STJ24751)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 IHC-P:1:50-1:200 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 101kDa Observed Mw: 101kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	4772
Gene Symbol	NFATC1
UniProt ID	NFAC1_HUMAN
Immunogen Region	20-300
Immunogen Sequence	VFGRGETLGPAPRAGGTMKS AEEEHYGYASSNVSPALPLP TAHSTLPAPCHNLQTSTPGI IPPADHPSGYGAALDGGPAG YFLSSGHTRPDGAPALESPI IEITSCGLYHNNNQFFHDV EVEDVLPSSKRSPSTATLSL PSLEAYRDPSCSPASSLSS RSCNSEASSYESNYSYPYAS PQTSPWQSPCVSPKTTDPEE GFPRGLGACTLLGSPRHSPS TSPRASVTEESWLGARSSR
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 20-300 of NFAT2 (NP_765978.1).

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunofluorescence analysis of Jurkat cells using NFAT2 Rabbit polyclonal antibody (STJ24751) at dilution of 1:100 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunohistochemistry analysis of NFAT2 in paraffin-embedded human placenta using NFAT2 Rabbit polyclonal antibody (STJ24751) at dilution of 1:200 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Western blot analysis of lysates from Mouse thymus, using NFAT2 Rabbit polyclonal antibody (STJ24751) at 1:600 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 180s.