



Anti-PDCD6IP antibody (1-180) (STJ24927)


GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/IF/ICC/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	Lot specific
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB:1:200-1:2000 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
	Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH7.3.
	Isotype	IgG
	Molecular Weight	Protein Mw: 96kDa Observed Mw: 105kDa
	Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	10015
	Gene Symbol	PDCD6IP
	UniProt ID	PDC6I_HUMAN
	Immunogen Region	1-180
	Immunogen Sequence	MATFISVQLKKTSEVDLAKP LVKFIQQTYPSGGEEQAQYC RAAEELSKLRRAAVGRPLDK HEGALETLLRYYDQICSIEP KFPFSENQICLTFTWKDAFD KGSLFGGSVKLALASLGYEK SCVLFNCAALASQIAAEQNL DNDEGLKIAAKHYQFASGAF LHIKETVLSALSREPTVDIS
	Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human ALIX/PDCD6IP (NP_037506.2).
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunofluorescence analysis of U-251MG cells using ALIX/ PDCD6IP antibody (STJ24927) at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of L929 cells using ALIX/ PDCD6IP antibody (STJ24927) at dilution of 1:100. Blue: DAPI for nuclear staining.



Western blot analysis of extracts of various cell lines, using ALIX/ PDCD6IP antibody (STJ24927) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 3min.