

Anti-PEBP1 antibody (1-187) (STJ24954)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/IP/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 IF/ICC:1:50-1:200 IP:0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 21kDa Observed Mw: 21kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	5037
Gene Symbol	PEBP1
UniProt ID	PEBP1_HUMAN
Immunogen Region	1-187
Immunogen Sequence	MPVDLSKWSGPLSLQEVD EQ PQHPLHVTYAGAAVDELGKV LTPTQVKNRPTSISWDGLDS GKLYTLVLTPDAPSRKDPK YREWHHFLVNMKGNDISSG TVLSDYVGSGPPKGTGLHRY VWLVEQDRPLKCDEPILSN RSGDHRGKFKVASFRKKYEL RAPVAGTCYQAEWDDYVPKL YEQLSGK
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-187 of human PEBP1 (NP_002558.1).

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Immunoprecipitation analysis of 150 Mu g extracts of MCF-7 cells using 3 Mu g PEBP1 antibody (STJ24954). Western blot was performed from the immunoprecipitate using PEBP1 antibody (STJ24954) at a dilution of 1:500.




Immunofluorescence analysis of PC-12 cells using PEBP1 Rabbit polyclonal antibody (STJ24954) at dilution of 1:25 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of NIH/3T3 cells using PEBP1 Rabbit polyclonal antibody (STJ24954) at dilution of 1:25 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of HeLa cells using PEBP1 Rabbit polyclonal antibody (STJ24954) at dilution of 1:25 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Western blot analysis of various lysates using PEBP1 Rabbit polyclonal antibody (STJ24954) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.