

Anti-PGAM1 antibody (1-254) (STJ24968)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:2000 IF/CC:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 29kDa Observed Mw: 29kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	5223
Gene Symbol	PGAM1
UniProt ID	PGAM1_HUMAN
Immunogen Region	1-254
Immunogen Sequence	MAAYKLVLRHGESAWNLEN RFSGWYDADLSPAGHEEAKR GGOALRDAGYEFDICFTSVQ KRAIRTLWTVLDAIDQMWLP VVRTWRLNERHYGGLTGLNK AETAAKHGAEQVKIWRRSYD VPPPPMEPDHPPFYNSISKDR RYADLTEDQLPSCESLKDTI ARALPFWNEEIVPQIKEGKR VLIAAHGNSLRGIVKHLEGL SEEAIMELNLP TGIPIVYEL DKNLKPDKPMQFLGDEETV
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-254 of human PGAM1 (NP_002620.1).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunofluorescence analysis of L929 cells using PGAM1 Rabbit polyclonal antibody (STJ24968) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Western blot analysis of extracts of various cell lines, using PGAM1 antibody (STJ24968) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.