

Anti-PGR antibody (1-250) (STJ24974)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/IHC-P/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	Lot specific
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB:1:500-1:1000 IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
	Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
	Isotype	IgG
	Molecular Weight	Protein Mw: 99kDa Observed Mw: 118kDa
	Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	5241
	Gene Symbol	PGR
	UniProt ID	PRGR_HUMAN
	Immunogen Region	1-250
	Immunogen Sequence	MTCLKAKGPRAPHVAGGPPS PEVGSPLLCPAAGPFPQSQ TSDTLPEVSAIPISLDGLLF PRPCGGQDPSEKTDQDQSL SDVEGAYSRAEATRGAGGSS SSPPEKDSGLLDVSLDTLLA PSGPGSQSPSPPACEVTSSW CLFGPELPEDPPAAPATQRV LSPLMSRSGCKVGDSSGTAA AHKVLPRGLSPARQLLLPAS ESPHWGAPVKPSPQAAAVE VEEEDGSESEESAGPLLKG
	Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human Progesterone Receptor (NP_000917.3).
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
   		