

Anti-TICAM1 antibody (180-712) (STJ25850)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 76kDa Observed Mw: 75kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	148022
Gene Symbol	TICAM1
UniProt ID	TICAM1_HUMAN
Immunogen Region	180-712
Immunogen Sequence	RPIDGVSDWSQGCSLRSTGS PASLASNLEISQSPTMPFLS LHRSPHGPSKLCDDPQASLV PEVPVGGCQEPPEMSWPPSG EIASPPPELPSSPPPGLEVA P
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 180-280 of human TRIF/TICAM1 (NP_891549.1).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of various lysates using TRIF/ TICAM1 Rabbit pAb (STJ25850) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution.

Lysates/proteins: 25 μ g per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit

Exposure time: 90s.