

Anti-TNXB antibody (444-673) (STJ25911)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 458kDa Observed Mw: 69kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7148
Gene Symbol	TNXB
UniProt ID	TNX_HUMAN
Immunogen Region	444-673
Immunogen Sequence	TSFTTGGLRIPFPRDCGEEM QNGAGASRTSTIFLNGNRER PLNVFCDMETDGGGWLVFQR RMDGQDFWRDWEDYAHGFG NISGEFWLGNEALHSLTQAG DYSMRVDLRAGDEAVFAQYD SFHVDSAAEYRHLHLEGYHG TAGDSMSYHSGSVFSARDRD PNSLLISCAVSYRGAWWYRN CHYANLNGLYGSTVDHQGVV WYHWKGFESVPFTEMKLRP RNFRSPAGGG
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 444-673 of human TNXB (NP_115859.2).

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Western blot analysis of various lysates using TNXB Rabbit polyclonal antibody (STJ25911) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.