

Anti-TRAF6 antibody (50-150) (STJ25955)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 60kDa Observed Mw: Refer to figures
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	7189
Gene Symbol	TRAF6
UniProt ID	TRAF6_HUMAN
Immunogen Region	50-150
Immunogen Sequence	MEEIQGYDVEFDPPLESKYE CPICLMALREAVQTPCGHRF CKACIIKSIRDAGHKCPVDN EILLENQLFPDNFAKREILS LMVKCPNEGCLHKMELRHLE D
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 50-150 of human TRAF6 (NP_004611.1).

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunofluorescence analysis of HeLa cells using TRAF6 rabbit polyclonal antibody (STJ25955) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Immunofluorescence analysis of NIH/3T3 cells using TRAF6 rabbit polyclonal antibody (STJ25955) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Immunofluorescence analysis of PC-12 cells using TRAF6 rabbit polyclonal antibody (STJ25955) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.