

Anti-WWP2 antibody (151-400) (STJ26122)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 99kDa Observed Mw: 105kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11060
Gene Symbol	WWP2
UniProt ID	WWP2_HUMAN
Immunogen Region	151-400
Immunogen Sequence	VPNGSALTDGSQLPSRDSSG TAVAPENRHQPPSTNCFGGR SRTHRHSGASARTTPATGEQ SPGARSRRHRQPVKNSGHSGL ANGTVNDEPTTATDPEEPSV VGVTSPPAAPLSVTPNPNTT SLPAPATPAEGEEPSTSGTQ QLPAAQAQAPDALPAGWEQRE LPNGRVYYVDHNTKTTTWER PLPPGWEKRTDPRGRFYVD HNTRTTTQWRPTAEYVRNVE QWQSQRNQLQGAMQHFSQR
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 151-400 of human WWP2 (NP_001257383.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of lysates from A-549 cells, using WWP2 Rabbit polyclonal antibody (STJ26122).