

## Anti-KIAA0391 antibody (314-583) (STJ26245)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Rat

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 67kDa Observed Mw: 67kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">9692</a>
Gene Symbol	<a href="#">PRORP</a>
UniProt ID	<a href="#">MRPP3_HUMAN</a>
Immunogen Region	314-583
Immunogen Sequence	YPGESFAHSIKTWFEVSPGK QWKGQFTTVRKSGQCSGCGK TIESIQLSPEEYECLKGKIM RDVIDGGDQYRKTTPELKR FENFIKSRPPFDVVIDGLNV AKMFPKVRSQLLNVSQ AKRNLRLVLGRKHMLRFS QWSRDEMEEVQKQASCFFAD DISEDPFLLYATLHSGNHC RFITRDLMRDHKACLPDAKT QRLFFKWQQGHQLAIVNRFP GSKLTFQRILSYDTVVQTT
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 314-583 of human KIAA0391 (NP_055487.2).

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of extracts of various cell lines, using KIAA0391 antibody (STJ26245) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 30s.