

Anti-EIF5A2 antibody (1-153) (STJ26921)

GENERAL INFORMATION

	Product Type	Primary antibodies
	Applications	WB/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:200-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 17kDa Observed Mw: 17kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	56648
Gene Symbol	EIF5A2
UniProt ID	IF5A2_HUMAN
Immunogen Region	1-153
Immunogen Sequence	MADEIDFTTG DAGASSTYPM QCSALRKNGFVVLKGRPCKI VEMSTSKTGKHGHAKVHLVG IDIFTGKKYEDICPSTHNMD VPNIKRNDYQLICIQDGYLS LLTETGEVREDLKLPEGELG KEIEGKYNAGEDVQVSMCA MSEEYAVAIPCK
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-153 of human EIF5A2 (NP_065123.1).

ADDITIONAL INFORMATION

Note: STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of extracts of various cell lines, using EIF5A2 antibody (STJ26921) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.