

Anti-WDFY4 antibody (1-160) (STJ26942)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 354kDa Observed Mw: 354kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57705
Gene Symbol	WDFY4
UniProt ID	WDFY4_HUMAN
Immunogen Region	1-160
Immunogen Sequence	MEAEADLSKAEDRNEDPGSKN EGQLAAVQPDVPHGGQSSSP TALWDMLEKRFLEYQQLTHK SPIERQKSLLSLLPLFLKAW EHSVGIICFPSLQRLAEDVS DQLAQQLKALVGKPAEQAR LAAGQLLWWKGDVDQDGYLL LKSVYVLTGTDSETLGRVAE
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human WDFY4 (NP_065996.1).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of lysates from HeLa cells, using WDFY4 Rabbit polyclonal antibody (STJ26942). Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST.