

Anti-TNFRSF10A antibody (269-468) (STJ28189)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 50kDa Observed Mw: 45kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8797
Gene Symbol	TNFRSF10A
UniProt ID	TR10A_HUMAN
Immunogen Region	269-468
Immunogen Sequence	GGDPKCMDRVCFWRLGLLRG PGAEDNAHNEILSNADSLST FVSEQQMESQEPADLTGVTV QSPGEAQCLLGPAEAEGSQR RRLLVPANGADPTETLMLFF DKFANIVPFDSDQLMRQLD LTKNEIDVVRAGTAGPGDAL YAMLMKWVNTGRNASIHTL LDALERMEERHAREKIQDLL VDSGKFIYLEDGTGSVLSLE
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 269-468 of human TNFRSF10A (NP_003835.3).

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Western blot analysis of extracts of DU145 cells, using TNFRSF10A antibody (STJ28189) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.